

EL

ESTHETIC LINE
implant



C-TECH
IMPLANT

TABLE OF CONTENTS

EL - ESTHETIC LINE IMPLANT

Implant characteristics	page 04
Dental implants	page 05
Closed tray impression transfer	page 06
Open tray impression transfer	page 07
Technical planning abutments	page 08
Titanium abutments	page 09
Temporary abutments	page 10
Screw-retained restorations	page 11
Closed tray technique	page 12
Screw retained restoration	page 13
Bar	page 14
O-Ball attachment system	page 15
Anchor abutment system	page 16
Deluxe Surgical Kit	page 17
Surgical Kit components	page 18
Drills & Bone Taps	page 19



All of the materials produced by C-TECH follow a validated procedure, which includes surface treatment and packing as well, in conformity with European and international directives EN ISO 13485:2003/AC:207 and 93/42/EEC relative to medical devices.



PRECISION DENTAL SOLUTIONS

C-Tech Implant is a dynamic company with aggressive growth, producing components and product lines primarily for dental implantology.

INTERNATIONAL PRESENCE

With production and management based in Italy, C-Tech Implant is active in all major world markets and is distributed in over 20 countries.

SCIENTIFIC RESEARCH, ADVANCED TECHNOLOGY, SIMPLIFICATION

C-Tech Implant differentiates itself with attention to research and the application of high technology to its products, all while maintaining a simplicity of insertion and ease of use.

C-Tech Implant incorporates the latest trends in implantology but provides very practical surgical and prosthetic solutions aimed at offering the practitioner and the patient optimal results.

HIGH QUALITY STANDARDS KEPT WITHIN REACH

C-Tech Implant products are made to the highest standards governing the manufacturing and management of European medical and dental components.

Up to date audits and certification assure that these standards are vigilantly maintained.

TRAINING & ADVICE

Dental professionals are assisted by the rich knowledge and experience of C-Tech Implant personnel and through C-Tech courses and training sessions.

During these courses the professional is able to learn the latest methods of implant placement and reconstruction.

MISSION STATEMENT

The goal of C-Tech Implant is to provide the highest level of quality for technologically advanced products at reasonable prices in order to allow the dental practitioner to find solutions for the broadest range of patients.



ESTHETIC LINE implant

BEVELLED SHOULDER

- Facilitates bone growth above the shoulder
- Long term implant stability
- Biological repartition of the forces in cortical bone

MICRO GROOVING

- Softens forces to the cortical bone during insertion
- Cortical bone maintenance

GRIT BLASTED AND ACID ETCHED SURFACE TOPOGRAPHY

- Best surface for osseo integration and bone to implant contact

AGGRESSIVE APICAL DESIGN

- Ideal for immediate implant placement
- Primary stability

ROUNDED APEX

- Protection of the sinus floor, nerve canal and other important anatomical structures during insertion

BONE LEVEL IMPLANT SUBCRESTAL SEATING

- Hinders exposure of the implant through bone resorption
- Ideal for the esthetic zone
- Long term esthetic stability

THREE DIFFERENT THREADING PROFILES

- Thread designs adapted to different bone structures that occur along the depth of the implant
- Enhanced surface
- Round but cutting apex design

DOUBLE LEAD THREAD

- Insertion rate of 1,5mm per rotation
- Primary stability
- Increased bone to implant contact
- Faster and even insertion while protecting bone structure

THREAD IN THREAD / GROOVE IN GROOVE

- Increased bone to implant contact

CONCAVE ESTHETIC CONCEPT

- Non surgical thickening of the peri-implantary tissue
- Facilitation of the papilla reconstruction-technique

PLATFORM SWITCHING

- Reduces bone loss
- Better representation of the biological width
- Long term esthetic stability

ONE CONNECTION FOR ALL 3 DIAMETERS

- Simplifies the system
- Reduces inventory
- Ease of use

COLD WELD SEAL


- Hinders bacterial infiltration and consequent bone loss

MORSE LOCKING CONICAL CONNECTION

- Elimination of micro-movements
- Elimination of screw loosening

INDEXING HEX


- Antirotational security



DENTAL IMPLANT



9 mm	EL-3509	EL-4309	EL-5109
11 mm	EL-3511	EL-4311	EL-5111
13 mm	EL-3513	EL-4313	EL-5113
15 mm	EL-3515	EL-4315	EL-5115



EL - SHORT IMPLANTS



7 mm	EL-4307	EL-5107
------	---------	---------

EL CEC PEEK HEALING ABUTMENTS Ø4

INCLUDES SCREW



EL-4502H



EL-4503H



EL-4504H*



EL-4506H*

EL CEC PEEK HEALING ABUTMENTS Ø5

INCLUDES SCREW



EL-5502H



EL-5503H



EL-5504H*



EL-5506H*



Please Note: the extractor screw (BL-6060 or BL-6061) is required to remove the healing abutment from the implant.

*Uses the long screw EL-5052HXL

CLOSED TRAY IMPRESSION TRANSFER



BL-4546
Peek Impression Cap



EL-4502P EL-5502P
EL-4503P EL-5503P
EL-4504P* EL-5504P*
EL-4506P* EL-5506P*

EL PEEK abutments together with the snap on BL-4546 cap function as impression transfers



EL-5502F EL-4502F
EL-5503F EL-4503F
EL-5504F* EL-4504F*
EL-5506F* EL-4506F*

EL Titanium Abutments together with the snap on BL-4546 cap function as impression transfers



BL-5143
Analog



BL-5052HX
Short Screw
EL-5052HXL
Long Screw

INTENDED USE

Closed tray impression technique.

CHARACTERISTICS

- Simple;
- Slender emergence profile to accommodate space limitations;
- No additional preparation (i.e. perforation) of tray required;
- High precision impression components give an exact replica of the intraoral situation;
- Clear-cut tactile response from the prosthetic connection verifies proper seating of components.

NOTE

Impression posts ensure optimal fit and precise impression taking for each patient.

* Uses long screw EL-5052HXL



STEP 1

Place the impression post accurately into the implant and hand-tighten the guide screw.

STEP 2

Push the impression cap at the top of the impression transfer.

STEP 3

Take the impression using an elastomeric impression material (polyvinyl siloxane or polyether rubber).

STEP 4

Use a standard impression tray.

STEP 5

Mount the impression transfer on the analog using the screw (ref. BL-5052HX - EL-5052HXL).

STEP 6

Reposition the impression transfer in the tray. Push the impression transfer until you feel the complete engagement firmly seated on the impression cap.

OPEN TRAY IMPRESSION TRANSFER



EL-4544
Open tray Impression
post includes the
BL-5050L guide screw



BL-5143
Analog



BL-5050L
Guide screw for open
tray impression post



BL-5050S
Short guide screw
for open tray
impression post

INTENDED USE

Open tray impression technique.

CHARACTERISTICS

- Simple;
- Slender emergence profile accommodates space limitations;
- Guide screw can be tightened either by hand or with the SCS screwdriver;
- High precision impression components give an exact replica of the intraoral situation;
- Clear-cut tactile response from the prosthetic connection verifies proper seating of components.

NOTE

Open tray impression procedure requires a custom-made tray with perforations. Impression posts are intended for single use only to ensure optimal fit and precise impression taking for each patient.



STEP 1

Place the impression post accurately into the implant and hand-tighten the guide screw.

STEP 2

Make perforations in the custom-made impression tray (light cured resin) according to the individual situation so that the positioning screw of the impression post sticks out.

STEP 3

Take the impression using an elastomeric impression material (polyvinyl siloxane or polyether rubber).


STEP 4

Reposition and fix the analog in the impression using the screw.

TECHNICAL PLANNING ABUTMENTS

Ø 4 STRAIGHT PLANNING ABUTMENTS

INCLUDES SCREW



INTENDED USE

Intra & extra-oral planning of prosthetic restoration.

CHARACTERISTICS


- Simple;
- Color-coded PLANNING abutments;
- Comprehensive PLANNING set containing all PLANNING abutments arranged clearly;
- Proper seating of PLANNING abutments verified through the clear-cut response from the prosthetic connection;
- PLANNING abutments fabricated of sterilizable polymer material.

NOTE

Be sure to clean and sterilize the planning abutments following intra-oral use. Do not sterilize the PLANNING abutment cassette.

Ø 5 STRAIGHT PLANNING ABUTMENTS

INCLUDES SCREW



STEP 1

Place the PLANNING abutment into the technical lab model situation in order to plan and choose the appropriate titanium abutment in cost effective manner.

STEP 2

Place the titanium abutment and hand-tighten the screw.

STEP 3

Prepare the titanium abutment, modify as required.

STEP 4



Fabricate the superstructure on the modified abutment using the standard modelling, casting and veneering methods.

STEP 5 - Cast the framework using the standard casting methods.

STEP 6 - Veneer the superstructure.

18° ANGLED PLANNING ABUTMENTS

INCLUDES SCREW



TITANIUM ABUTMENTS

EL CEC TITANIUM Ø4 ABUTMENTS

INCLUDES SCREW



EL-4502F


EL-4503F

EL-4504F*

EL-4506F*

EL CEC TITANIUM Ø5 ABUTMENTS

INCLUDES SCREW



EL-5502F


EL-5503F

EL-5504F*

EL-5506F*

EL CEC 18° ANGLED ABUTMENTS

INCLUDES SCREW



EL-4518/2

EL-4518/3

EL SCAN ABUTMENT



EL-SCAN
Compatible with EXOCAD,
3SHAPE and
DENTALWINGS

TITANIUM/ZIRCONIUM BASE



BL-6046

Complete set includes titanium base and prosthetic screw

TITANIUM CASTABLE ABUTMENT



BL-6045

Complete set includes titanium base, casting cylinder and screw

TITANIUM CEREC® BASE



BL-6047

Complete set includes titanium cerec® base and prosthetic screw

INTERNAL PROSTHETIC SCREWS



BL-5052HX
Standard
prosthetic screw



EL-5052HXL
Long prosthetic
screw for H4 - H6

INTENDED USE

Cement-retained restorations.

TIGHTENING:


with torque ratchet 30
N=Ncm

CHARACTERISTICS

- Simple;
- Less grinding necessary due to prepared mucosa margins;
- Adaptation to natural soft tissue contour due to prepared mucosa margins in different heights (H2, H3, H4, H6);
- Reliable;
- Tapered connection (pure cone). Abutment and implant are linked so as to form a one-piece unit;
- Extractor system allows easy abutment removal from the implant or the analog.

NOTE

The cement margin must not be more than 2 mm below the mucosa. Use a new basal screw for the final insertion of the abutment.



ABUTMENT EXTRACTOR SCREW

As the ABUTMENT EXTRACTOR SCREW is driven in, it will push the abutment out of the analog or implant.



CT-E7003
Finger/Ratchet
adapter for latch
drivers



BL-6060
Prosthetic extractor




BL-6061
Latch driver
prosthetic extractor

*Uses the long screw EL-5052HXL

TEMPORARY ABUTMENTS


EL CEC PEEK Ø 4 ABUTMENTS

INCLUDES SCREW AND IMPRESSION CAP



EL CEC PEEK Ø 5 ABUTMENTS

INCLUDES SCREW AND IMPRESSION CAP



*Uses the long screw EL-5052HXL

NOTE

Together with the BL-4546 Impression Cap, the EL temporary PEEK abutments can be used as closed tray impression transfers.



BL-4546
Peek Impression
Cap

NOTE

Together with the EL-4543 Temp-Cap, the EL temporary PEEK abutments can be used to temporarily stabilize a prosthesis.



EL-4543
Peek Temp-Cap

ZIRCONIUM ABUTMENTS

All zirconium abutment packagings include screw.

ZIRCONIUM STRAIGHT ABUTMENTS



ZIRCONIUM 15° ANGLED ABUTMENTS




ZIRCONIUM 25° ANGLED ABUTMENTS




Please Note: prosthetic screw for zirconium abutment reference is EL-5052HXL

SCREW-RETAINED RESTORATIONS



STRAIGHT ABUTMENTS



17° ANGLED ABUTMENTS



30° ANGLED ABUTMENTS




COMPLETE SET INCLUDES:



1. Multi-unit angled abutment
2. Prosthetics Screw



CLOSED TRAY TECHNIQUE



STEP 1
Remove the healing abutments.



STEP 2
Screw the straight abutment into the implant using the torque ratchet (30 Ncm) and the Multi-unit Driver.

SURGICAL PROCEDURE



STEP 1
Screw each closed tray transfer onto the protruding abutments.



STEP 2
Take the impression using an elastomeric impression material (polyvinyl siloxan or polyether rubber).



STEP 3
Remove the closed tray transfer from the abutment.

STEP 4
Screw onto the abutments the healing cap screws so as to keep the soft tissue in place until the final prosthesis is completed.

LABORATORY PROCEDURE



STEP 1
Screw the closed tray transfer onto the analog.

STEP 2
Reposition the transfer into the previously taken impression material being sure that the transfers are properly seated.




STEP 3
Master model.

SCREW RETAINED RESTORATION

①



②



STEP 1

Fabricate the stone model including analogs and gingival mask.

STEP 2

Place and screw the castable abutments onto the protruding multi-unit analogs.

STEP 3

Shorten the cylinders down to the height of the occlusal plane.

③



④



STEP 4

Remove the gingiva modeling material to permit easy access for submucosal contouring and verification of component seating. Wax-up the bridge framework to appropriate dimensions. The layer of wax must have sufficient thickness to avoid the wrong coefficient of thermal expansion and a negative effect on porcelain firing.

STEP 5

Prepare the wax-up for investing and casting procedures.

STEP 6

Attach the resulting framework to the models and create final prosthesis.

STEP 7

Passively fit the resulting prosthesis onto the abutments.

⑤




⑥



BAR



BL-5146
Multi-unit Analog



0220BB
OT-Bar (2 pcs.)



027CRR
Clip pink: soft (4 pcs.)




027CRG
Clip yellow: medium
(4 pcs.)



BL-5647
Castable
Abutment
includes screw




BL-6051
Bridge screw



STEP 1
Place the castable Multi-unit
abutments on the analogs and tighten
the Multi-unit internal screws.

STEP 2
Make height adaptations according to
the individual situation.




STEP 3
Use a residue-free burn-out plastic to
fix the bar segments to the castable
abutments.

STEP 4
The yellow clips (027CRG) are fixed
into the prosthesis.

O-BALL ATTACHMENT SYSTEM

O-ball
Abutment Driver



INTENDED USE

Removable dentures retained by implants in the mandible and maxilla.

CHARACTERISTICS

- Simple;
- The clinical process for the ball attachment is quick and easy;
- Functional;
- The O-ring attachment is designed to virtually eliminate wear on the Ball Abutment and minimize the need for maintenance;
- 3 different gingival heights;
- 3 different O-ring resistances offering optimal retention for every individual situation.

RELIABLE

Dual retention for optimal abutment-denture connection. Excellent long-term performance due to wear resistant components.

STEP 1

Screw the spherical abutment into the implant using the torque ratchet (30 Ncm) and the driver (ref. BL-0600).


STEP 2

Rebase the overdenture according to standard procedure.

STEP 3

Use a laboratory burr to relieve the denture base in the indicated areas.

METAL HOUSING



COMPLETE SET INCLUDES:



1. O-Ring (Ref. MC-3005 , MC-3005B)
1 piece
2. Metal Housing (Ref. MCH-2)
3. O-Ball Abutment (Ref. 5644, 5642, 5641)

ANCHOR ABUTMENT SYSTEM

CAPS WITH METAL HOUSING



141CAE
2 Stainless steel housings



140CEV
4 Retentive caps
violet "strong" (2.7kg)



140CET
4 Retentive caps
white "standard" (1.8kg)




140CER
4 Retentive caps
pink "soft" (1.2kg)



140CEG
4 Retentive caps
yellow "extra-soft" (0.6kg)

ANCHOR ABUTMENT SYSTEM COMPLETE SET



130BL1



130BL2



130BL3



130BL4



130BL5



- 1 Anchor abutment
(Ref. 130BL1, 130BL2, 130BL3, 130BL4, 130BL5)
1 Stainless steel housings (Ref. 141CAE)
1 Retentive caps - violet "strong" (Ref. 140CEV)
1 Retentive caps - white "standard" (Ref. 140CET)
1 Retentive caps - pink "soft" (Ref. 140CER)
1 Retentive caps - yellow "extra-soft" (Ref. 140CEG)

LABORATORY ACCESSORIES



140CEN
4 Processing Caps - black



144MTE
2 Impression Coping



144AE
2 Laboratory Analog



044CAIN
Pull-off Impression Coping

SURGICAL INSTRUMENTS



124ICP
1 Blue plastic
"multiuse"
insertion tool



185IAC
1 Metal insertion tool for
caps



191ECS
1 Metal extractor tool for
caps



774CQ
1 OT-Equator square
screw driver for implant
abutment (square 1,25mm)



760CE
1 Square driver
connector for torque

DELUXE SURGICAL KIT

DELUXE SURGICAL KIT INCLUDES:

Locator drill CT-2020	Depth Gauge (3.5 mm) CT-E9007
2.0 external irrigation drill (Ø 2.0) CT-1720E	Depth Gauge (4.3 mm) CT-E9008
3.5 external irrigation drill (Ø 3.0) CT-1735E	Depth Gauge (5.1 mm) CT-E9010
4.3 external irrigation drill (Ø 3.6) CT-1743E	Paralleling Pins, qty. 2 (1.6 mm & 2.0 mm) CT-9000
5.1 external irrigation drill (Ø 4.6) CT-1751E	Set metal stopper (L.9/11/13/15) CT-Stop07/08/09/10
Hard Bone Drill 3.5 mm EL-1735N	Implant Latch Driver BL-E9040
Hard Bone Drill 4.3 mm EL-1743N	Implant Ratchet Driver Short BL-E7001
Hard Bone Drill 5.1 mm EL-1751N	Implant Ratchet Driver Long BL-E7001L
Drill Extender CT-2000	Extractor BL-6060
Manual Hex Driver Short CT-9025S	Implant ratchet driver (short) ND-E7001
1.25mm Hextool Torque Wrench Attachments CT-8051	Implant ratchet driver (LONG) ND-E7001L
1.25mm Hextool Torque Wrench Attachments (Long) CT-8052	3.0 external irrigation drill ND-1726E
Torque Wrench (50 Ncm) CT-8010	

REFERENCES:

EL-Sur Kit.01:
Deluxe Surgical Kit for EL Line



SURGICAL KIT COMPONENTS

SURGICAL FACILITATORS

COMBINED DEPTH PARALLELING PINS



CT-PIN1
1.9 MM
2.5 MM

PARALLELING PINS



CT-9000
1.6 mm
2.0mm

DEPTH GAUGES



CT-E9007
3.5mm
ø est.



CT-E9008
4.3mm
ø est.



CT-E9010
5.1mm
ø est.



CT-2000
Drill Extender

PROSTHETICS DRIVER



CT-9025S
Hex drivers



CT-9025
Hex drivers



CT-9019
Handpiece Latch



CT-8051
Torque wrench
attachments



CT-8052
Torque wrench
attachments



BL-6060
Prosthetic extractor



BL-6061
Latch driver
prosthetic extractor

DRILL DEPTH STOPS



CT-STOP06
Stop L.6



CT-STOP02
Stop L.7



CT-STOP01
Stop L.8



CT-STOP07
Stop L.9



CT-STOP03
Stop L.10



CT-STOP08
Stop L.11



CT-STOP12
Stop L.12



CT-STOP09
Stop L.13



CT-STOP14
Stop L.14



CT-STOP10
Stop L.15

IMPLANT DRIVERS

(WITH RETENTION)



BL-E7001
Implant
ratchet driver



BL-E7001L
Implant
ratchet driver



BL-E9040
Implant latch

FINGER ADAPTORS



CT-E7002
Finger adaptor
for ratchet drivers



CT-E7003
Finger/ratchet adaptor
for latch drivers


TORQUE WRENCH 50Ncm

CT-8010



TORQUE WRENCH PEEK

CT-8010PEEK



DRILLS & BONE TAPS



SURGICAL PROCEDURE EL Ø 3.5 IMPLANT



CT-2020
Locator Drill



CT-1720E
2.0mm Irr. Drill
(Outer Ø 2.0mm)



CT-1735E
3.5mm Irr. Drill
(Outer Ø 3.0mm)



EL-1735N
3.5 mm
hard bone drill




SURGICAL PROCEDURE EL Ø 4.3 IMPLANT



CT-2020
Locator Drill



CT-1720E
2.0mm Irr. Drill
(Outer Ø 2.0mm)



CT-1735E
3.5mm Irr. Drill
(Outer Ø 3.0mm)



CT-1743E
4.3mm Irr. Drill
(Outer Ø 3.6mm)



EL-1743N
4.3 mm
hard bone drill



SURGICAL PROCEDURE EL Ø 5.1 IMPLANT



CT-2020
Locator Drill



CT-1720E
2.0mm Irr. Drill
(Outer Ø 2.0mm)



CT-1735E
3.5mm Irr. Drill
(Outer Ø 3.0mm)



CT-1743E
4.3mm Irr. Drill
(Outer Ø 3.6mm)



CT-1751E
5.1mm Irr. Drill
(Outer Ø 4.6mm)



EL-1751N
5.1 mm
hard bone drill

***IMPORTANT:**

0.8 mm must be added to the length of the drill considering for the angled cutting tip



Via Santa Margherita al Colle n. 18 - 40136, Bologna - ITALY
Tel. +39 051 6661817 - Fax +39 051 6667071
www.c-tech-implant.com - info@c-tech-implant.com

Go to
www.c-tech-implant.com



Follow us on facebook
[f/CTechImplant](https://www.facebook.com/CTechImplant)

